


# Biopotential electrodes

- Bioelectric signals
- Biopotential electrode models
- Electrode types

Biomedical Electronics and Instrumentation  
Department of Electronic Engineering

# Bioelectric activity


- ◆ Electrical signals produced by various bodily activities are used in monitoring / diagnosis
  - Measured signals are generally compared with normal reference signals
  - Variations in signals can provide accurate diagnosis
- ◆ Examples include
  - Electrocardiography (ECG) – activity of different sections of the heart
  - Electroencephalography (EEG) – electrical activity of the brain
  - Electromyography (EMG) – electrical activity of the muscles
- ◆ Signals are characterised by
  - Amplitude
  - Frequency / time duration
  - Waveshape



**Characteristics of  
ECG waveforms**

# Biopotential signals


- ◆ Electrical activity is explained by differences in ion concentrations within the body (sodium,  $\text{Na}^+$ ; chloride,  $\text{Cl}^-$ ; potassium,  $\text{K}^+$ )
  - A potential difference (voltage) occurs between 2 points with different ionic concentrations
  - Cell membranes at rest are more permeable to some ions (e.g.  $\text{K}^+$ ,  $\text{Cl}^-$ ) than others (e.g.  $\text{Na}^+$ )
    - $\text{Na}^+$  ions are pumped out of the cells, while  $\text{K}^+$  ions are pumped in
    - Due to a difference in rates of pumping ( $\text{Na}^+ > \text{K}^+$ ), a difference in +'ive ion concentration results
    - A -'ive electrical potential is established between the inside and outside of the cell; this is estimated to be approx.  $-70 \text{ mV}$ .
    - Cell is said to be *polarised*



## ◆ Action potential


- When a cell is electrically stimulated (through the central nervous system), the permeability of the cell membrane changes
  - $\text{Na}^+$  ions rush into the cell, and  $\text{K}^+$  ions rush out
  - Again, due to a difference in rates of flow, the ion concentration changes (more +’ive ions inside cell than outside)
  - Cellular potential becomes +’ive : 20 – 40 mV. This lasts for a few milliseconds
  - Cell is said to be *depolarised*
- After the stimulus, the permeability of the cell membrane returns to its original value, and the rest potential is restored
  - Due to unequal pumping rates of  $\text{Na}^+$  and  $\text{K}^+$  ions
  - Time taken for restoration is called the refractory period
  - Cell is said to be *repolarised* during this time
- The resulting variation in cellular potential with time is known as the action potential

# Typical action potential signal



- ♦ ECG, EEG etc signals are the resultant of several action potentials produced by a combination of different cells

# Action potentials → measured waveforms



- ◆ Propagation of action potential through different body tissues produces final waveform recorded by electrodes

# Biopotential sensors / electrodes

- ◆ Biopotential electrodes transduce *ionic* conduction to *electronic* conduction so that biopotential signals can be viewed and/or stored
- ◆ They generally consist of metal contacts packaged so that they can be easily attached to the skin or other body tissues



Metal plate electrodes



Suction electrode





Needle electrodes

- Different types include surface macroelectrodes, indwelling macroelectrodes & microelectrodes
- Skin and other body tissues act as electrolytic solutions

# Half-cell potential

- ◆ The performance of an electrode may be described in terms of a metal – electrolyte interface




- Electrode discharges some metallic ions into electrolytic solution
  - Increase in # free electrons in electrode and in # +ive cations in solution
- OR ions in solution combine with metallic electrodes
  - Decrease in # free electrons in electrode and in # +ive cations in solution
- As a result, a charge gradient builds up between the electrode and electrolyte and this in turn creates a potential difference; i.e. the electrode potential or half-cell potential,  $V_h$ .



# Electrode double layer

- ◆ For both mechanisms, two parallel layers of oppositely charged ions are produced; i.e. the electrode double layer :
  - E.g. when metal ions recombine with the electrode :




- The excess of negative anions is replaced with positive cations in the case of metal ions discharging into solution, and  $V_h$  is then  $< 0$ .

# Polarisation

- ◆ A perfectly polarised electrode is one in which there is no net transfer of charge across the metal / electrolyte interface.
  - Either metal ions are discharged OR electrolytic ions recombine
  - Only displacement current flows
  - Interface impedance is represented as a capacitor
- ◆ A perfectly non-polarised electrode is one in which there is no impedance to charge transfer between the metal and electrolyte
  - Both mechanisms of charge transfer occur
  - There is no double layer build up
  - Interface impedance is represented as a resistor
- ◆ In practice, both mechanisms occur to some extent and interface impedance therefore has components of resistance *and* capacitance



# Electrode circuit model



- ◆ Biopotential is measured as the difference in readings between two electrodes
  - Ideally, the half-cell potential is cancelled out
- ◆ Model assists in amplifier design (later)
  - Source resistance vs. amplifier input impedance
  - Frequency dependance

# Typical electrode materials

reduction reaction	$E^{\circ}$ (V)
$Al^{3+} + 3e^{-} \rightarrow Al$	-1.662
$Zn^{2+} + 2e^{-} \rightarrow Zn$	-0.762
$Cr^{3+} + 3e^{-} \rightarrow Cr$	-0.744
$Fe^{2+} + 2e^{-} \rightarrow Fe$	-0.447
$Cd^{2+} + 2e^{-} \rightarrow Cd$	-0.403
$Ni^{2+} + 2e^{-} \rightarrow Ni$	-0.257
$Pb^{2+} + 2e^{-} \rightarrow Pb$	-0.126
$2H^{+} + 2e^{-} \rightarrow H_2$	0.000 ← Standard Hydrogen electrode
$AgCl + e^{-} \rightarrow Ag + Cl^{-}$	+0.222
$Hg_2Cl_2 + 2e^{-} \rightarrow 2Hg + 2Cl^{-}$	+0.268
$Cu^{2+} + 2e^{-} \rightarrow Cu$	+0.342
$Cu^{+} + e^{-} \rightarrow Cu$	+0.521
$Ag^{+} + e^{-} \rightarrow Ag$	+0.780
$Au^{3+} + 3e^{-} \rightarrow Au$	+1.498
$Au^{+} + e^{-} \rightarrow Au$	+1.692


*Note: Ag-AgCl has low junction potential & it is also very stable -> hence used in ECG electrodes!*

# Problem of half-cell potential

- ◆ Typical half-cell potential values (at DC) are much larger than AC measurement signals
  - $V_e : 0.1 - 0.5 \text{ V}$  vs.  $V_{\text{ECG}} : 1 - 2 \text{ mV}$  &  $V_{\text{EEG}} : \sim 50 \mu\text{V}$
  - Even if identical electrodes are used and a differential signal is measured, there will be differences in  $V_e$  values due to variation in placement positions
- ◆ Solutions include :
  - Use of a differential DC amplifier – however, small remaining differences will be amplified by the same level as low level signals to be detected.
  - Include a DC offset circuit to subtract the half-cell potential – this approach is limited by the time-varying nature of the offset due to motion between the skin and electrode.
  - Include a high pass filter to remove the DC offset – this provides the best results for measurement signals with frequency much higher than that of offset variations

# Surface electrodes


- ◆ These are placed in contact with the skin of the subject



- ◆ Problem of motion artifact
  - Electrode movement produces unwanted signals due to change in electrode – electrolyte interface


# Floating electrodes

- ◆ Conductive paste reduces effect of electrode slippage and resulting motion artifact :



# Needle electrodes

- ◆ Mostly used for contacting with internal body tissues :




- (a) Insulated needle electrode .
- (b) Coaxial needle electrode .
- (c) Bipolar coaxial electrode .
- (d) Fine -wire electrode connected to hypodermic needle, before being inserted .
- (e) Cross -sectional view of skin and muscle, showing coiled fine -wire electrode in place .



# Microelectrodes

- ◆ Used to measure biopotential signals at the cellular level :

## Metal Supported Microelectrodes




(a) Metal inside glass


(b) Glass inside metal

- ◆ Due to small dimensions ( $\mu\text{m}$ ), impedance levels are high
  - Amplifier needs very high input impedance

# Microtechnology / MEMS



(a) Beam -lead multiple electrode



(b) Multielectrode silicon probe